



Diagnostic Slides

In this leaflet you will find a selection of our slides with printed reaction fields for the medical diagnostic (e.g.: immunofluorescent, histology or other microbiological tests). 26 x 76 x 1.0 mm, according to ISO 8037/1, high quality soda-lime float glass, ground edges 45°, with or without twin frosted end. PTFE or epoxy, standard black or other colours (white, red, green, blue etc.). Cleanly washed and polished, ready to use, packed by 50 pieces in fibre-free cardboard boxes. Glass surface either wettable or with adhesive coating. Custom slides according to your specification are also available.











Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	60 / Ø 2 mm	K 001
	20 mm one side	36 / Ø 2 mm	K 002
	20 mm one side	15 / Ø 6 mm	K 003
	20 mm one side	15 / Ø 6 mm	K 004
	20 mm one side	12 / Ø 6 mm	K 005
	20 mm one side	12 / Ø 6 mm	K 006
	20 mm one side	12 / Ø 6 mm	K 106
	20 mm one side	12 / Ø 5 mm	K 007











Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	12 / Ø 5 mm	K 008
	20 mm one side	12 / Ø 5 mm	K 108
	20 mm one side	12 / Ø 4 mm	K 009
	20 mm one side	12 / Ø 4 mm	K 010
	20 mm one side	10 / Ø 8 mm	K 011
	20 mm one side	10 / Ø 7 mm	K 012
	20 mm one side	10 / Ø 6 mm	K 013
	20 mm one side	10 / Ø 6 mm	K 014











Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	10 / Ø 5 mm	K 015
	20 mm one side	10 / Ø 5 mm	K 016
	20 mm one side	8 / Ø 9 mm	K 017
	20 mm one side	8 / Ø 8 mm	K 018
	20 mm one side	8 / Ø 6 mm	K 019
	20 mm one side	8 / Ø 6 mm	K 020
	20 mm one side	6 / Ø 10 mm	K 021
	20 mm one side	6 / Ø 10 mm	K 022











Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	6 / Ø 8 mm	K 023
	20 mm one side	6 / Ø 8 mm	K 024
	20 mm one side	5 / Ø 9 mm	K 025
	20 mm one side	5 / Ø 8 mm	K 026
	20 mm one side	5 / Ø 4 mm	K 027
	20 mm one side	5 / 21,9 x 10 mm	K 028
	20 mm one side	4 / Ø 12 mm	K 029
	20 mm one side	4 / Ø 10 mm	K 030











Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	4 / Ø 9 mm	K 031
	20 mm one side	4 / Ø 8 mm	K 032
	20 mm one side	3 / Ø 15 mm	K 033
	20 mm one side	3 / Ø 11 mm	K 034
	19 mm one side	3 / 18 x 25 mm	K 035
	20 mm one side	2 / 19 x 19 mm	K 036
	10 mm two sides	3 / Ø 10 mm	K 037
	20 mm one side	6 / Ø 7 mm	K 038











Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	3 / Ø 8 mm	K 039
	without	8 / Ø 6 mm	K 040
	20 mm one side	2 / Ø 15 mm	K 041
	20 mm one side	2 / Ø 11 mm	K 042
	20 mm one side	2 / Ø 8 mm	K 043
	20 mm one side	2 / Ø 6 mm	K 044
	20 mm one side	1 / Ø 20 mm	K 045
	20 mm one side	1 / Ø 15 mm	K 046











Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	1 / Ø 12 mm	K 047
	20 mm one side	1 / Ø 10 mm	K 048
	20 mm one side	1 / Ø 8 mm	K 049
	20 mm one side	1 / Ø 6 mm	K 050
	without	80 / 2 x 2 mm	K 051
	without	4 / Ø 10 mm	K 052
	without	4 / 9 x 9 mm	K 053
	without	5 / Ø 6 mm	K 054











Image	Frosted End	Fields / Measures	K-Nr.
	without	30 / Ø 4 mm	K 055
	20 mm one side	24 / Ø 4 mm	K 155
	without	24 / Ø 5 mm	K 056
	20 mm one side	18 / Ø 5 mm	K 156
	without	16 / Ø 6 mm	K 057
	without	16 / Ø 3 mm	K 058
	without	12 / Ø 6 mm	K 059
	without	8 / Ø 12 mm	K 060











Image	Frosted End	Fields / Measures	K-Nr.
	without	6 / Ø 10 mm	K 061
	without	4 / 10 x 10 mm	K 062
	without	4 / Ø 10 mm	K 063
	without	3 / Ø 15 mm	K 064
	20 mm one side	2 / Ø 11 mm	K 065
	20 mm one side	4 / Ø 11 mm	K 066
	20 mm one side	4 / Ø 6 mm	K 067
	without	16 / Ø 5 mm	K 068











Image	Frosted End	Fields / Measures	K-Nr.
	33 mm one side	1 / 41 x 24 mm	K 069
	20 mm one side	4 / Ø 8 mm	K 070
	20 mm one side	1 / Ø 11 mm	K 071
	20 mm one side	2 / Ø 6 mm	K 072
	20 mm one side	3 / Ø 11 mm	K 073
	20 mm one side	3 / Ø 10 mm	K 074
	10 mm one side	12 / Ø 8 mm	K 075
	10 mm one side	12 / Ø 8 mm	K 076




















Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	9 / 5 x 20 mm	K 077
	20 mm one side	1 / Ø 15 mm	K 078
	15 mm one side	10 / Ø 4 mm	K 079
	20 mm one side	16 / Ø 5 mm	K 080
	20 mm one side	8 / Ø 6 mm	K 081
	20 mm one side	30 / Ø 2 mm	K 082
	20 mm one side	6 / Ø 4 mm	K 083
	without	1 / Ø 12 mm	K 085



Image	Frosted End	Fields / Measures	K-Nr.
	20 mm one side	18 / Ø 4 mm	K 086
	10 mm one side	10 / Ø 8 mm	K 087
	without	5 / Ø 8 mm	K 088
	10 mm one side	10 / Ø 8 mm	K 089
	20 mm one side	1 / Ø 13 mm	K 090
	20 mm one side	12 / Ø 5 mm	K 095
	without	8 / Ø 11 mm	K 096
	20 mm one side	14 / Ø 4 mm	K 097





Contact:

W. Knittel Glasbearbeitungs GmbH
Varrentrappstraße 5
38114 Braunschweig
Germany

Fon: +49 (0)531 / 59 000-0
Fax: +49 (0)531 / 59 000-36
Email: info@knittelglass.com
Web: www.knittelglass.com

