

029



**Obstetrics**  
**Geburtshilfe**  
**Obstetricia**  
**Obstétrique**  
**Ostetricia**



½

**DE LEE**  
 29-104-20  
 20 cm, 8"



½

**COLLIN**  
 29-108-35  
 35 cm, 13¾"



½

**MARTIN**  
 29-110-34  
 34 cm, 13½"



½

**COLLYER**  
 29-114-30  
 30 cm, 12"



½

**BREISKY**  
 29-116-33  
 33 cm, 13"



½

**THOLE**  
 29-118-34  
 34 cm, 13¾"



¼

**SIMPSON**  
 29-120-23  
 23 cm, 9"



¼

**SIMPSON-BRAUN**  
 29-124-30 30 cm, 12"  
 29-124-36 36 cm, 14"




¼

**SIMPSON-LUIKART**  
 29-130-36  
 36 cm, 14"



¼

**DE LEE**  
 29-136-30 30 cm, 12"  
 29-136-36 36 cm, 14"  
 29-136-39 39 cm, 15½"



¼

**McLEAN**  
 29-140-36  
 36 cm, 14"



¼

**NAEGELE**  
 29-142-36 36 cm, 14"  
 29-142-40 40 cm, 16"



¼

**KIELLAND**  
 29-146-41  
 41 cm, 16"



¼

**KIELLAND-LUIKART**  
 29-148-40  
 40 cm, 15½"



$\frac{1}{4}$   
**WRIGLEY**  
29-150-28  
28 cm, 11"



$\frac{1}{4}$   
**PIPER**  
29-154-44  
44 cm, 17 $\frac{1}{4}$ "



$\frac{1}{4}$   
**TARNIER**  
29-158-40  
40 cm, 15 $\frac{3}{4}$ "




$\frac{1}{4}$   
**BARNES-NEVILLE**  
29-160-43  
43,5 cm, 17"



$\frac{1}{4}$   
**SELLHEIM**  
29-170-28  
28 cm, 11"

**Obstetric lever**  
Geburtsöffel  
Separador para la salida de cabeza  
Levier obstétrique  
Leva ostetrica



$\frac{1}{4}$   
**MURLESS**  
29-175-30  
30 cm, 12"

**Head extractor for use in lower segment in caesarean operation**

Kopfnabel für unteres Segment bei Kaiserschnitt

Elevador para cabeza usándose en le segmento inferior del útero en casos de operaciones cesáreas

Levier à tête utilisé en cas de césarienne

Elevatore per testa in interventi con taglio cesareo



$\frac{1}{2}$   
**BILL**  
29-177-00

**Traction Handel to be used with all forceps with standard type handle**

Zeghebel, passend für alle Geburtszangen mit gehorntem Standard-Griff

Tractor abecudado para todas las forceps con mango standard

Tracteurs pour tous les forceps avec manches ordinaires

Manico a trazoine adatto a tutte le forcipi con manico standard




**GREEN-ARMYTAGE**  
 29-180-21  
 21 cm, 8 1/4"

**Caesarean section haemostatic forceps**  
 Uterusklemmzange zur Blutstillung  
 Pinza-clamp uterina para hemostasis  
 Pince-clamp utérine pour hémostase  
 Pinza emostatiche per tagli cesarei



**WILLET**  
 29-184-19  
 19 cm, 7 1/2"

**Scalp flap forceps**  
 Kopfschwartenzange  
 Pinza para pericráneo  
 Pince pour le cuir chevelu  
 Pinza per il cuoio cappelluto



29-190-01



29-190-02



29-190-03



29-191-01



29-191-02



29-191-03

**WINTER**  
 29-190-01  
 29-191-03  
 28 cm, 11"

**Without retchet**  
 Ohne Sperre  
 Sin cremallera  
 Sans crémaillère  
 Senza cremagliera



**SAENGER**  
 29-199-27  
 27,5 cm, 10 3/4"

**With retchet**  
 Mit Sperre  
 Con cremallera  
 Avec crémaillère  
 Con cremagliera


**SAENGER**  
 29-198-27  
 27 cm, 10 1/2"

**Placenta and ovum scoops**

Plazentar- und abortuslöffel  
Cucharillas para aborto y placenta  
Curettes à avortement et placenta  
Cucchiai per placenta e ovaie

**Placenta and ovum curettes**

Plazentar- und Abortusküretten  
Curetetas para aborto y placenta  
Curettes à avortement et placenta  
Curette per placenta e ovaie



¼  
29-200-01 24 mm  
29-202-01 24 mm



¼  
29-200-02 26 mm  
29-202-02 26 mm



¼  
29-200-03 28 mm  
29-202-03 28 mm

½  
**BUMM**  
29-200-01 29-202-01  
29-200-03 29-202-03  
31 cm, 12¼" 31 cm, 12¼"

**Blunt, Malleable**  
Stumpf, biegsam  
romas, maleable  
mousses, malléable  
Smusse, malleabile

**Sharp, malleable**  
Scharf, biegsam  
Cortantes, maleable  
Trachentes, malléables  
acute, malleabile



¼  
29-206-35  
35 mm



¼  
29-206-40  
40 mm

½  
**RECAMIER**  
29-206-35  
29-206-40  
31 cm, 12¼"


**Placenta and ovum curettes,  
blunt / rigid**

Plazentar- und Abortusküretten,  
stumpf / fest

Curetetas para aborto y placenta,  
romas / rígidas

Curettes à avortement et placenta,  
mousse / rigide

Curette per placenta e ovaie,  
smusse / rigide



½  
**CUZZI**  
29-210-01  
29-210-03  
30 cm, 12"

**Placenta and ovum scoops**  
Plazentar- und abortuslöffel  
Cucharillas para aborto y placenta  
Curettes à avortement et placenta  
Cucchiai per placenta e ovaie



¼  
29-210-01  
14,5 mm



¼  
29-210-02  
18 mm



¼  
29-210-03  
21 mm



1/2

**SPELLIE**  
29-216-25  
25,5 cm, 10"


**Perforator**  
Perforatorium  
Perfora-crâneos  
Perce-crâne  
Spezzacrano



1/2

**BRAUN**  
29-220-42  
42 cm, 16 1/2"

**Cranioclast**  
Kranioklast  
Craneoclasto  
Cranioclaste  
Cranioclasto



1/2

**BRAUN**  
29-224-32  
32,5 cm, 12 3/4"

**Decapitating hook**  
Dekapitationshaken  
Gancho para decapitación  
Crochet pour décapitation  
Uncino per detroncazione



1/2

**COLLIN**  
29-230-09  
9 cm, 3 1/2"

**Umbilical cord clamp**  
Nabelschnurklemme  
Clamp umbilical  
Clamp pour le cordon ombilical  
Clamp ombelicale



1/2

**HASSELLTINE**  
29-232-06  
6 cm, 2 3/4"

**Umbilical cord clamps, in packages of 12 pieces**

Nabelschnurklemmen, in packungen zu 12 Stück

Clamps umbilicales, en paquetes de 12 piezas

Clamps pour le cordon ombilical, en paquet de 12

Pinze ombelicali, pacchi da 12 pz.



1/2

**IOWA**  
29-236-14  
14 cm, 15 1/2"

**Trumpet needle guide**  
Nadelführung  
Aguja de guía  
Aiguille trompette guide  
Cuida ago a tromba

